

Monitoring intracellular protein interactions using NanoLuc® Binary Technology (NanoBiT™)

Brock Binkowski, Christopher Eggers, Braeden Butler, Marie Schwinn, Michael Slater, Thomas Machleidt, Mei Cong, Keith Wood & Frank Fan


Promega Corporation, 2800 Woods Hollow Rd, Madison, WI 53711-5399



1. Introduction


Protein:protein interactions (PPIs) are essential to the cellular signal transduction pathways that contribute to cancer. Although numerous approaches exist to monitor PPIs *in vitro*, methods for intracellular detection have been more limited. We developed NanoLuc® Binary Technology (NanoBiT), a two-subunit system based on NanoLuc® luciferase that can be applied to the intracellular detection of PPIs. Large BiT (LgBiT; 18 kDa) and Small BiT (SmBiT; 11 amino acid peptide) subunits are expressed as fusions to proteins of interest, where PPI facilitates subunit complementation to give a bright, luminescent enzyme. Unlike related approaches where an enzyme or protein is simply split, LgBiT was independently optimized for structural stability and SmBiT was selected from a peptide library specifically for the PPI application. The result is a subunit pair that weakly associates ($K_D = 190 \mu\text{M}$) yet still maintains 30% of the activity of full-length NanoLuc at saturation. In contrast to many split systems, the LgBiT:SmBiT interaction is reversible, allowing the detection of rapidly dissociating proteins. PPI dynamics can be followed in real-time in living cells using the Nano-Glo® Live Cell Reagent, a non-lytic detection reagent containing the cell-permeable furimazine substrate. Advantages over split systems include better sensitivity, reversibility, fusion to a peptide or a small, structurally stable protein domain, real-time measurements using a non-lytic assay format, and subunits with reduced affinity for self-association.

2. NanoBiT overview




- LgBiT and SmBiT are fused to proteins A & B
- A:B interaction facilitates LgBiT:SmBiT interaction, generating a bright luminescent enzyme
- LgBiT:SmBiT with low affinity ($K_D = 190 \mu\text{M}$), limiting non-specific association and reducing assay background
- LgBiT:SmBiT interaction is reversible ($k_{on} = 500 \text{ M}^{-1}\text{sec}^{-1}$; $k_{off} = 0.2 \text{ sec}^{-1}$)
- LgBiT evolved for increased structural stability making it a better fusion partner
- Non-lytic assay format allows real-time measurements of protein interaction dynamics for 1-2 hrs

3. Real time kinetics of PPI interaction dynamics




- Transient expression of SmBiT-PRKACA & LgBiT-PRKAR2A in HEK293
- Modulators of intracellular cAMP added sequentially at indicated time points
- Inverse correlation for NanoBiT vs. GloSensor cAMP 22F (cAMP biosensor)
- NanoBiT can monitor reversible PPIs in real time


4. GPCR interaction with β -arrestin-2



5. RAF dimerization using stable expression




6. Nuclear hormone receptor dimerization




- Transient expression of LgBiT-AR & AR-SmBiT in HEK293 cells
- R1881 agonist added at time zero, which induces AR dimerization
- CRC data plotted at 30 minutes
- Validated for GR homodimerization as well
- NanoBiT can monitor nuclear hormone receptor dimerization in real time

7. Agonist-induced RTK dimerization



7. Quantifying agonist binding via dimerization



9. Summary

NanoBiT is extremely bright

- Fusion partners can be expressed at very low levels, minimizing potential artifacts
- 100-1,000 fold brighter than split firefly luciferase

NanoBiT components are small & stable

- LgBiT, 18 kDa; SmBiT, 11 amino acids
- LgBiT evolved for increased structural stability, providing a stable fusion partner

NanoBiT is reversible

- Monitor both protein association and dissociation events in real time

NanoBiT offers experimental flexibility

- Monitor protein interaction dynamics at a single time point or continuously for 1-2 hours
- Room temperature or 37 °C measurements
- Validated in 96-, 384- & 1536-well formats

Check <http://www.promega.com/nanobit> for new NanoBiT PPI expression vectors

For more information on NanoBiT, please see *ACS Chemical Biology*, 11(2), p. 400-408